

Contents

Section	Page
<u>Airway Care</u>	4 - 16
<u>Bio Safety Supplies</u>	17 - 23
<u>Branded & Custom</u>	24 - 27
<u>Crash and Jump bags</u>	28 - 36
<u>Diagnostics and Monitoring</u>	37 - 64
<u>Disinfectants</u>	65 - 68
<u>Disposables</u>	69 - 95
<u>Emergency Medical Services and Safety Equipment</u>	96 - 104
<u>First Aid Equipment</u>	105 - 114
<u>First Aid Kits</u>	115 - 127
<u>Immobilisation Devices</u>	128 - 135
<u>Respiratory Care</u>	136 - 147
<u>Stretchers and Lifting Devices</u>	148 - 161
<u>Vascular Access & Care</u>	162 - 166
<u>Be Safe Paramedical Contacts</u>	167
<u>Order Form</u>	168

Click on link to be taken to corresponding page

AIRWAY CARE

Version 6-2013

Four Blade Laryngoscope Set

Product Code: LARYN.COMP

- Manufactured from high quality Japanese non-magnetic stainless steel
- Single piece construction ensures blade integrity
- Smooth satin finish makes cleaning and sterilisation optimal
- Set comes with size 0 and 1 Miller Blade and size 3 and 4 Macintosh Blade
- Handle design promotes comfortable grip
- Durable soft carry case included



Four Blade Laryngoscope Set – Fibre Optic

Product Code: LARYN.FO.COMP



- Manufactured from high quality non-magnetic Japanese stainless steel
- Single piece construction ensures blade integrity
- Smooth satin finish makes cleaning and sterilisation optimal
- Set comes with size 0 and 1 Miller Blade and size 3 and 4 Macintosh Blade
- Removable and replaceable fibre optic bundles
- Handle design promotes comfortable grip
- Durable hard carry case included

Five Blade Laryngoscope Set – Fibre Optic

Product Code: LARYN.COMP

- Lightweight and compact
- Manufactured from high quality non-magnetic Japanese stainless steel
- Smooth satin finish makes cleaning and sterilisation optimal
- Removable light guides for easy cleaning
- 3 spring ball-bearings provides secure fit and stability
- Set includes size 0 Miller blade and size 1, 2, 3 and 4 Macintosh blades
- Adult and paediatric handle included



Six Blade Laryngoscope Set – Fibre Optic

Product Code: LARYN.FO.PB.COMP



- Blades manufactured from high quality plastic
- Set comes with size 2, 3 and 4 Macintosh blade and size 1, 2 and 3 Miller blade
- Blades are single use
- Handle design makes for a comfortable grip and excellent hygiene
- Durable hard carry case included

Oropharyngeal Airway Box Set – Complete

Product Code: AB.COMP

- Set of eight oropharyngeal airways contained in transparent plastic container
- Airway sizes include: 000, 00, 0, 1, 2, 3, 4 and 5.
- Manufactured from high quality plastic
- All edges on airways are rounded to reduce the risk of oral trauma
- Airways colour coded for easy identification
- Airways are for single and individual use



Nasopharyngeal Airways

Product Code: NASO.AIR.



- Excellent option for managing the airway in the semi-conscious patient
- Manufactured from medical grade PVC
- Large flange minimises the chance of tube slippage
- Soft flexible texture reduces the chances of nasal trauma during insertion
- Available in singles and in a full range of sizes
- Individually packaged

Endotracheal Tubes – Sizes 2.5 – 9.5 (Cuffed)

Product Code: ET.TC

- Manufactured from non-toxic, clear thermo sensitive siliconised PVC material
- Cuff with integral one-way-valve
- All tubes fitted with 15mm connector compatible with standard ventilation devices
- Full length radio-opaque line to identify tube placement
- Tube marking in increments of 1cm
- 2cm mark assists with placement past vocal cords



Endotracheal Tubes – Sizes 2 – 5.0mm (uncuffed)

Product Code: ET.TUC



- Manufactured from non-toxic, clear thermo sensitive siliconised PVC material
- All tubes fitted with 15mm connector compatible with standard ventilation devices
- Full length radio-opaque line to identify tube placement
- Tube marking in increments of 1cm
- 2cm mark assists with placement past vocal cords

Laryngeal Mask Airway – Sizes 1 – 5

Product Code: LARYN.A.MASK

- Manufactured from high quality PVC
- Clear tube to allow for early detection of fluids
- Latex free
- Highly visible non-removable blue ink printing for easy measurement
- High quality 15mm connector
- Classic cuff shape
- Touch sensitive pilot balloon for accurate cuff pressure measurement



I-Gel Supraglottic Airway

Product Code: AIR.I-GEL



- Made from a medical grade thermoplastic elastomer
- designed to create a non-inflatable, anatomical seal of the pharyngeal, laryngeal and perilaryngeal structures whilst avoiding compression trauma
- Available in sizes 1- 5
- Individually packaged
- Single use

Endotracheal Tube Stylette (ETT Introducer)

Product Code: INTRO

- Provides rigidity to endotracheal tubes to facilitate endotracheal intubation
- Manufactured from metal covered medical plastic
- Available in various sizes
- Sterile
- Single use



Endotracheal Tube Introducer (Bougie)

Product Code: INTRO.BOUG.I.



- Used to facilitate difficult endotracheal intubation
- Ultra smooth surface for low friction
- Markings per cm
- Angled coude tip for ease of intubation
- Supplied sterile, individually packaged in easy to open sachets
- Available in 15, 10 and 6Fr sizes to cover full range of patients

CritiTube Check (Oesophageal Detection Device (ODD))

Product Code: CRITI.TC.I

- The CritiTube Check is used to rapidly assess endotracheal tube placement following intubation
- Cost-effective option to disposable CO₂ detectors
- 15mm adapter
- Illuminates in dark
- Available in yellow



Endotracheal Tube Holder – Adult

Product Code: ET.TH.T.A



- Used to secure endotracheal tube following intubation
- Extremely easy to apply
- Provides excellent ETT stability
- Adult size covers 5.0 to 10mm range of endotracheal tubes
- Built in bite block removes need for use of an oropharyngeal airway
- Foam padding on inside of holder promotes patient comfort
- Very cost-effective to use

CritiPack® Surgical Airway Pack

Product Code: CRITI.P.SA.I

- All equipment required to effectively perform a surgical airway via endotracheal tube
- Pre-packed in transparent packaging
- Sterile
- Fully disposable
- High quality products
- Checklist skills guide for procedure included
- Excellent focus on hygiene and patient safety



CritiJet Jet Insufflation Device

Product Code: CRITI.J



- Cricothyroidotomy performed in the “can’t intubate, can’t ventilate” situation can be life-saving
- The CritiJet developed and manufactured by Be Safe Paramedical is used to provide jet insufflation after performing needle cricothyroidotomy
- Extremely easy to use and far less cumbersome than expensive jet insufflation devices
- Luer lock system provides secure connection to IV catheter and extra length tubing allows the provider to position himself/ herself behind the patient during operation
- Provides excellent pressures

Emergency Surgical Airway Instrument Set

Product Code: EMERG.N.TRACH

- All instruments manufactured from high quality medical stainless steel
- Set includes; scalpel, scalpel blade, scalpel handle, scissors, bb tissue forceps, toothed tissue forceps, non-toothed tissue probe scissors
- Instruments contained in zip up bag for secure storage



Pocket Mask Resuscitator

Product Code: RES.M.PKT



- Manufactured from high quality transparent PVC
- Latex free
- Pre-inflated cushion for quick effective seal
- One-way-valve incorporates a filter for added infection protection
- Standard fittings allows device to be used with a bag-valve resuscitator
- Suitable for adult and paediatric use
- Oxygen inlet allows for high-flow oxygen to be connected to mask
- Nose marking for reduced risk of misplacement

Suction Catheter with Control Port - Yankhauer

Product Code: SUC.YANK.WC

- Manufactured from high quality clear PVC
- Non-toxic and non-pyrogenic
- Latex free
- Catheter end designed for atraumatic suctioning
- Finger control port for controlled suctioning application
- Individually packed



Suction Catheter with Control Port –Soft Tip

Product Code: SUC.C.WC



- Manufactured from high quality clear PVC
- Non-toxic and non-pyrogenic
- Latex free
- Available in sizes; 5FG, 6FG, 8FG, 10FG, 12FG, 14FG, 16FG and 18FG
- Ends of catheters are designed for atraumatic suctioning
- Thumb control port for controlled suctioning application
- Colour-coded port for easy size identification
- Individually packed

Suction Unit – Zeiner Rescue (Battery operated)

Product Code: SUC.RESCUE

- Size: 240 x 125 x 216 mm
- Mass: 2,8 Kg (including battery)
- Current: 3,5 A
- Adjustable vacuum up to 610 mmHg
- Max flow: 25 l/min
- Battery 3300 mAh (Ni-MH battery cells)
- Supplied with bag
- Warranty: 2 years
- Made in Italy



Suction Unit – Zeiner Rescue Plus (Battery operated)

Product Code: SUC.RESCUE.PL



- Size: 240 x 125 x 216 mm
- Mass: 2,8 Kg (including battery)
- Current: 3,5 A
- Run from A/C or DC supply
- Max flow: 25 l/min
- Battery 3300 mAh (Ni-MH battery cells)
- Ambulance mounting bracket
- Direct electrical and mechanical connection to bracket
- Can be detached from bracket and used as a portable unit
- Warranty: 2 years
- Made in Italy

Suction Unit – H003 Portable Battery Operated

Product Code: SUC.H003

- Lightweight and fully portable
- Oil free lubrication pump
- Pressure adjustment control
- Removable 1000ml fluid reservoir
- Low noise
- Voltage: AC220V 50Hz
- Max negative pressure value: $\geq 0.070\text{MPa}$
- Pumping rate: 18L/min
- Power: 90VA
- Runs off A/C or internal battery



Suction Unit – Hand Held Manual

Product Code: SUC.HNDHELD



- Manufactured from high quality plastic
- Compact design easily allows for storage in jumpbag
- Provides 500mmHg vacuum
- 300ml disposable container
- Smooth surface promotes hygienic use
- Comfortable trigger grip for prolonged use

[Back to Contents Page](#)

Bio Safety Supplies

Bio Suit

Product Code: BIO.SUIT

- Manufactured from 3 layer plastic foil for maximum protection and durability
- Protects the user against direct contact with contaminants
- Elasticised wrist and ankle straps
- FFP3 respiratory mask included
- Green Nitrile Rubber gloves included
- Uvex 9301 protective goggles included
- Suit can be decontaminated and reused



Bio Suit

Product Code: BIOBAG.EBV.30



- Patient bed within unit is a vacuum mattress and fully adjustable to the contours of patient's body
- The cockpit of the Bio-Bag is attached to the vacuum mattress by means of two zips. This allows for complete detachment for storage and decontamination
- Canopy manufactured using high frequency welding technology, which promotes strength, durability and quality
- Respiratory filter built into patient canopy
- Built in battery is rechargeable

Complete Biological Protective System

Product Code: BIO.COMP

- Complete system for use in biological disaster situations
- Includes decontamination unit, bio suit and bio bag
- All manufactured from high quality materials
- Contact us for a brochure



The Mobile Workplace

Product Code: MOB.WORK



- Tents are manufactured from polyamide textile, which is extremely strong and durable. The material is also easy to disinfect
- Tent floors are manufactured from double coated polyvinylchloride textile. An additional plastic cover can be fitted to floors, which promotes easy disinfecting and high level hygiene
- Inflatable tent frame is double chambered with an over pressurisation valve
- Tents can be connected together in any combination by means of a universal connecting sleeve

Sharps Containers

Product Code: SHARPS.950ML

- Manufactured from quality plastic
- Used for the safe disposal of sharps such as needles, scalpel blades and suture cutters
- Well marked to ensure safe handling and storage
- Clip down lids
- Available in 950ml, 5Lt, 7.6Lt, 10Lt, 20Lt and 25Lt sizes



Sharps Container – Jump Bag/Desktop Model (200ml)

Product Code: SHARPS.200ML



- Manufactured from quality plastic
- Used for the safe disposal of sharps such as needles, scalpel blades and suture cutters
- Well marked to ensure safe handling and storage
- Compact size and design makes this an excellent sharps container to keep in your jump bag or on your desk to reduce the chances of a needle stick injury

Surgical Gloves – High Risk Blue

Product Code: GLOV.HR.L

- Manufactured from extra thick latex which is chlorinated and pH balanced to provide superior chemical resistance and protection
- Available in S, M, L and XL sizes
- Boxes of 100



Surgical Gloves – Latex

Product Code: GLOV.100.L



- Single-use examination glove, ambidextrous
- 100% latex
- Lightly powdered with corn starch
- Hypoallergenic.
- High tactile sensitivity

Surgical Gloves – Nitrile

Product Code: GLOV.N.L.PR

- Made of 100% Nitrile, a synthetic rubber which contains no Latex
- Puncture resistant
- Powder-free
- Suitable for latex-sensitive users
- Provides comfortable secure fit
- Available in S, M, L and XL
- Available as pairs or boxes of 100



Surgical Mask – N95

Product Code: MASK.DUST



- 3-ply design respirator produced from high filterable and electrostatic material
- Latex-free head straps
- Adjustable nose clip
- Low respiratory resistance with high filtration efficiency
- Hygienic, environmentally friendly and non-irritating to skin

Surgical Mask – Ear loop or tie

Product Code: MASK.SUR

- 3 ply latex-free
- No glass fibres
- Hypoallergenic
- Low resistance to breathing
- High filtration capacity
- Excellent fit
- Available with ear loop or tie



[Back to Contents Page](#)

Branded & Custom

Examples of our Custom Developed and Branded First Aid Kits for various organisations.

Contact us for all your custom designed corporate First Aid Kits

The Outdoor Warehouse Range



The Red Cross Range



The Sportmans Warehouse Range



The St. John Range



The Telkom Range



[Back to Contents Page](#)

Crash & Jump Bags

CritiBag - 1st Responder Jump Bag

Product Code: BAG.1ST

- Manufactured from high quality nylon
- Double stitched central attachment points
- Heavy duty zips used in all internal access points
- Comes standard with 7 internal removable storage pouches
- Internal pouches are colour-coded
- Transparent covers on internal pouches make for easy identification of contents
- Bag can be carried as a shoulder sling or as a hand pack



CritiBag – Medipack Jump Bag

Product Code: BAG.MP.I



- Manufactured from high quality nylon
- Double stitched central attachment points
- Heavy duty zips used in all internal access points
- Comes standard with 6 internal removable storage pouches
- Internal pouches are colour-coded with transparent covers
- Can be carried as a back pack or hand pack

CritiBag - Trauma

Product Code: BAG.T

- Manufactured from high quality nylon
- Double stitched central attachment points
- Heavy duty zips used in all internal access points
- Comes standard with 4 internal removable storage pouches
- Internal pouches are colour-coded with transparent covers
- Can be carried as a back pack or hand pack



CritiBag – Intermediate Life Support

Product Code: BAG.ILS.I



- Manufactured from high quality nylon
- Double stitched central attachment points
- Heavy duty zips used in all internal access points
- Comes standard with 11 internal removable storage pouches
- Internal pouches are colour-coded with transparent covers
- Removable Intravenous Access bag included
- Can be carried as a back pack or hand pack

CritiBag – Advanced Life Support

Product Code: BAG.ALS

- Manufactured from high quality nylon
- Double stitched central attachment points
- Heavy duty zips used in all internal access points
- Comes standard with 15 internal removable storage pouches
- Internal pouches are colour-coded with transparent covers
- Removable drug bag, airway bag and diagnostic bag included
- Can be carried as a back pack or hand pack



CritiBag – Paramedic Drug Bag (Large)

Product Code: BAG.D.P.L



- Manufactured from high quality nylon
- Insulated covers
- Double stitched zip
- 16 small pouches with zips in red section
- 4 large pouches with zips in blue section
- Blue section lockable for safe keeping of controlled medications
- Document pouch and pen holder inside cover

Mini Grabber Bag

Product Code: BAG.GM

- The baby brother of the legendary Grabber Bag from Be Safe Paramedical
- Manufactured from high quality nylon
- Heavy duty zip at main access point
- 5 internal removable storage pouches
- Storage pouches have transparent covers
- Double handles make for comfortable carrying



Grabber Bag

Product Code: BAG.G



- The Legendary Grabber Bag designed and manufactured by Be Safe Paramedical is one of the most widely used smaller bags around
- Manufactured from high quality nylon
- Heavy duty zip at main access point
- 6 internal removable storage pouches
- Storage pouches have transparent covers
- Double handles make for comfortable carrying
- More packing volume than the Mini Grabber

Doctor's Bag

Product Code: BAG.DOCTORS

- Excellent option for medical practitioners requiring more space for additional items of equipment
- Manufactured from high quality nylon
- Convenient access through central zip and flip open lid
- Internal compartments make for convenient and safe storage of items
- 2 small removable storage bags included
- Can accommodate a good volume of equipment



Contents for demo purposes only and are excluded. A fully stocked back option is available...

First Aid Bag 107

Product Code: BAG.107



- Excellent little bag with good volume to accommodate basic first aid equipment
- Manufactured from high quality nylon
- An excellent option to be branded with your company logo and packaged according to your needs
- Dimensions: 27cm x 30cm x 3cm

First Aid Bag 109
Product Code: BAG.109

- More volume than the 107 First Aid Bag and therefore more packing space
- Manufactured from high quality nylon
- Once again, an excellent option to be branded with your company logo and packaged according to your needs
- Dimensions: 27cm x 33cm x 9.5cm



First Aid Bag 140
Product Code: BAG.140



- Manufactured from high quality nylon
- As with all our bags, can be custom branded and packed according to your needs
Small enough to place in your travel bag with some basic supplies when on that next business or holiday trip
- Dimensions: 25cm x 16cm x 2.5cm

First Aid Bag 310

Product Code: BAG.140

- More packing volume than the 140 bag and suitable for a slightly extended packing list of items
- Manufactured from high quality nylon
- Dimensions: 28cm x 27cm x 2.5cm
- Can be custom branded and packed according to your needs



First Aid Bag - Moonbag

Product Code: MOONBAG



- Excellent for use on the sports field for those extras you need quick access to
- Compact design ensures a comfortable fit and won't get in the way when kneeling or bending
- Manufactured from high quality nylon
- Can be custom branded and packed according to your needs
- Dimensions: 10cm x 19cm x 30cm

Be Safe Paramedical IFAK Bag – Leg Rig

Product Code: BAG.IFAK.LR.BL

- Manufactured from high quality nylon
- Internal storage strapping provides for a broad range of equipment
- No branding or reflective tape, thus making bag excellent for the tactical and military environments where space limitations are crucial
- This rig attaches to a belt via a clip system and then is strapped to the upper thigh by means of a buckle and clip
- Available in black or camo
- Available as a bag only or fully equipped kit



Be Safe Paramedical IFAK Bag – Belt Rig

Product Code: BAG.IFAK.BR.BL



- Manufactured from high quality nylon
- Pull-out bag with nylon bungy cord
- High volume compact design
- Attaches to belt or back pack via straps and clip
- Available in black or camo
- Available as a bag only or fully equipped kit

[Back to Contents Page](#)

Diagnostics & Monitoring

Alcohol Breath Tester – DA-3000

Product Code: ABT.DA-3000

- One button operation and superb accuracy makes this alcohol breath tester an excellent choice for the entry level user
- Semiconductor sensor
- 3 digit LCD display
- Runs off 2 AAA batteries
- Warm up time: 20 seconds
- Response time: 5 seconds



Alcohol Breath Tester – DA-7000

Product Code: ABT.DA-7000



- Advanced high grade semiconductor
- 4 digit display with backlight
- Units of measurement can be adjusted
- Warm up time: 20 seconds
- Response time: 5 seconds
- Displays total number of tests performed
- Calibration every 6 months or 200 tests
- Soft carry case included

Alcohol Breath Tester – DA8000

Product Code: ABT.DA8000

- Fuel Cell Sensor
- Sensor cleaning mechanism activates during power-up
- CE and RoHS certified
- Sensor lasts at least 500 samples before replacement required. (when used according to guidelines set out in instruction manual)
- $\pm 0.005\% \text{BAC}$ at $0.100\% \text{BAC}$ When a person's actual breath alcohol concentration is 0.100% , the DA8000 will always produce a result between 0.105% to 0.095%
- Adjustable Sampling Time improves accuracy by ensuring the breathalyzer captures the entire breath sample
- Measurements in mg/l, %BAC and g/l
- Dieters and Diabetics users are known to produce a high level of acetone in their breath which the fuel cell sensor can determine which would otherwise be wrongly detected as alcohol
- Shock-absorbing case helps to protect the sensitive sensor from accidental shocks which may adversely affect the accuracy of the breathalyzer



Alcohol Breath Tester – DA-9000

Product Code: ABT.DA9000



- This is the newest and most cutting edge addition to our range of alcohol breath testers
- Enhanced accuracy by means of an advanced fuel cell sensor
- Printer enabled, with printer available as an optional extra
- Warm up time: 18 seconds
- Response time: 10 seconds
- Over 500 readings on 2 AA batteries
- Units of measurement can be adjusted
- Calibration every 500 tests or 6-12 months

Alcohol Breath Tester – DA-9000 Hand Held Printer

Product Code: ABT.DA9000.P

- Wireless mobile thermal printer for the DA-9000 alcohol breath tester



Sphygmomanometer - Aneroid Deluxe

Product Code: BPM.ANER



- Comfortable and durable adult cuff with Velcro fastener
- Artery positioning label
- Loop for gauge attachment
- Standard air release valve
- Leatherette zip case included

Sphygmomanometer - Aneroid with stethoscope

Product Code: BPM.AN.ST

- Calibrated nylon cuff with 'D' ring
- Double tube adult sized latex bladder
- Standard latex bulb
- Air release valve with spring
- Standard end valve
- Single head stethoscope



Sphygmomanometer – Palm Type

Product Code: BPM.AN.SH



- Palm type for single handed operation
- 50mm chrome plated gauge
- Standard adult blue vinyl cuff
- Standard latex bulb and valve
- Zippered leatherette case

Sphygmomanometer – Multi Cuff Set

Product Code: BPM.P.M5C



- The best solution for covering your full range of patients in one compact kit
- Set contains: thigh cuff, obese cuff, standard adult cuff, paediatric cuff and a newborn cuff
- Standard air release valve
- All cuffs made from blue durable nylon
- Rubber gauge protector
- Supplied in a quality compact, marked blue nylon bag

Sphygmomanometer - Aneroid Wall/Desk

Product Code: BPM.A.W

- Two variants in one design (desk or wall)
- Square-shaped with a 14cm diameter dial
- Front and rear adjustable for easy interpretation
- Standard adult cotton cuff
- Double tube adult sized latex bladder
- Coiled PVC tubing



Sphygmomanometer – Cuffs

Product Code: BP.C.



- Neonatal, Infant, Paediatric, Adult, Obese and Thigh sizes available in both single and double tube options and in various designs
- Compatible with most manual and electronic vital signs monitors
- All manufactured from high quality materials
- Reliable and cost effective option to replace your damaged or used BP cuff

Sphygmomanometer – Digital Wrist

Product Code: BPM.W.DIG.BP600W

- Liquid crystal display
- Displays pulse rate
- Takes 4 AA batteries
- BP Range 30-280mmHg
- Pulse range 40-200 Bpm
- Displays year, date and time



Sphygmomanometer – Digital Arm

Product Code: BPM.DIG.IMP.BPCB0A2H



- Liquid crystal display
- Displays pulse rate
- Oscillometric measuring method
- Low battery indicator
- Auto power off: 180 seconds

Reanibex Series 200 Automated External Defibrillator

Product Code: DEF.REAN.S200.AED

- Biphasic technology - Truncated Exponential Waveform adapted to patient's impedance
- Energy delivered: up to 200 joules
- Lithium battery with more than 300 discharge capacity
- Audible and visual prompts to guide users through resuscitation sequence
- Rugged and durable for use in all types of environments
- Portable and lightweight for easy transportation
- Optional PC software to store and review ECG and audio data
- Automatic daily self-test ensures the unit is always ready
- Configurable CPR protocol



Accessories

- Carrying case
- Adult disposable Cable-Electrodes for defibrillation
- Adult Pre-connected disposable cable-electrodes
- Non rechargeable battery
- Reanibex Data Management Software



Options

- Voice recording and storage of audio data
- ECG signal display



Reanibex Series 300 Automated External Defibrillator

Product Code: DEF.REAN.S300.AED



- Biphasic AED or AED with Manual override option
- Output energy: Up to 360 Joules
- Rugged and compact for easy use in any environment
- Large colour display 5,7"
- User friendly interface
- Universal AED algorithm for adult and paediatric patients
- ECG monitoring (Optional)
- Synchronised Cardioversion (on manual option)
- Automatic daily self-test ensures that it is always ready to use
- CPR guidance according to international protocols
- USB connection
- Two battery options: non-rechargeable and rechargeable
- Auto-train function – allows unit to be used as an AED trainer, therefore no need for separate trainer



Accessories

- Carrying case
- Ambulance bracket (EN 1789 certified)
- 3-lead patient cable (optional for unit with manual override)
- Adult multifunction electrodes
- Paediatric multifunction electrodes
- Rechargeable battery/ Non rechargeable battery
- ECG signal display through external leads

Reanibex Series 500 Automated External Defibrillator

Product Code: DEF.REAN.S500.AED

- Manual and Automated External Biphasic Defibrillator up to 360 Joules
- Lightweight with monitoring capabilities
- Rugged and compact for use in any environment
- Large colour display 5,7"
- User friendly ergonomic design
- Universal AED algorithm for adult and paediatric patients
- Synchronised Cardioversion
- 3-lead patient cable (optional)
- Monitor (IEC 60601-2-27)
- Auto-train function – unit can be used as an AED trainer without need for separate unit
- Automatic daily self-test ensures the device is always ready to use
- CPR guidance according to international protocols
- Non-invasive pacing (optional)
- Pulse Oximetry with Masimo SET technology (optional)
- USB connection
- Non- rechargeable and rechargeable battery options
- Resuscitation guidelines: Factory set guideline 2010 (ERC/AHA)

Accessories

- Carrying case
- Ambulance bracket (EN 1789 certified)
- CE certified
- Adult multifunction electrodes
- Paediatric multifunction electrodes



Reanibex Series 700 Cardiac Monitor/ Defibrillator

Product Code: DEF.REAN.S700.WP



- Biphasic technology
- Truncated Exponential Waveform adapted to patient's impedance
- Simple and intuitive
- Portable and lightweight
- Control of the defibrillator and the printer from the paddles
- Integrated adult/ paediatric paddles
- Operates from the mains (AC) Vehicle Battery (DC) and from its internal battery
- Autotest when switched on and during operation
- Internal and external flash memory card (on unit with AED option)
- CE and EN 1789 certified



Accessories

- 4-lead patient cable
- 5-lead patient cable (standard)
- 10-lead patient cable
- Internal paddles
- Multifunction electrodes
- Carrying case

Options

- Pulse oximetry
- External Pacing
- "ECG Visor" software
- "EL" electroluminescent screen



Reanibex Series 800 Cardiac Monitor/ Defibrillator

Product Code: DEF.REAN.S800M.WP

- Manual and Automated Biphasic Defibrillate up to 360 Joules
- Lightweight with monitoring and diagnostic capabilities
- Large colour display (8.4")
- High performance with a focus on ease of use and portability
- Integrated adult/ paediatric paddles (standard)
- Modular system - easily upgradeable to include various options – no need to for technician, modules are simply slotted in
- 12-lead ECG data transmission from the field to the Emergency Service or receiving emergency unit
- Bluetooth transmission (optional)
- Synchronized Cardioversion
- 3, 5 (standard) and 10- lead patient cable
- Masimo Rainbow SET technology that monitors SpO2, Carboxyhemoglobin and Methemoglobin (optional)
- Oridion capnography - Microstream technology (optional)
- Non-invasive Blood Pressure from SunTech Medical (optional)
- Operates from the mains (AC), from a vehicle battery (DC) and from its dedicated internal battery
- Various options and accessories available
- CE and EN1789 certified
- Read Dot Design Award for 2013



Automated External Defibrillator – AED 7000

Product Code: DEFIB.AED7000



- Great as a public access AED
- 3-step defibrillation process
- 2-button operation
- Extensive verbal and visual prompts during operation
- Continuous event recording which can either be printed or downloaded to a PC
- Weekly self-test top ensure operational readiness
- Biphasic energy output
- Not suitable for paediatric patients
- Lock-out protection to prevent inadvertent defibrillation
- Accessories available



UT100C EtCO₂/ SPO₂ Monitor

Product Code: VSM.UT100C.H.CONF:A

- The UT100C with its reliable technology and unbelievable pricing makes it the most affordable portable Capnograph on the SA market
- Colour LCD screen with screen rotation
- Auto power off function for saving battery life
- Three selectable operational modes
- Long battery life – up to 20 hours on 4 AA Alkaline batteries or Ni-MH rechargeable batteries
- Powerful data storage capacity
- Patient information management
- Data can be transferred to PC for storage, review and printing
- Audio and visual alarm capability
- CE Certified



Technical Specifications:

- Dimensions: 13.5 x 7.5 x 2.8cm
- Weight: 258g with batteries
- Battery Type: Alkaline or Ni-MH rechargeable
- >20 hours with <6 hours charge time
- Available in various configurations, pulse oximetry, capnography or combination capnography/ pulse oximetry

Accessories

- Protective cover
- Range of pulse oximeter probes
- Ni-MH batteries and external charger
- Adult and paediatric airway adapters



UT100 Hand Held Pulse Oximeter

Product Code: PO.HH.UT100



- 2.8 inch display with screen rotation
- Reliable measurement of SpO2 and Pulse Rate
- Trend data storage time up to 120 hours
- Visual and audible alarms
- Three working modes
- Four display modes
- Data can be transferred to PC for analysis review and printing
- Various probe options available
- Only 258 grams with cells
- CE certified

NT1A Hand Held Pulse Oximeter

Product Code: PO.HH.UT100

- Designed for spot-check or continuous monitoring
- Audible/visible Alarms and adjustable alarm limits
- Digital SpO2 technology
- Accurate Measurement During Low Perfusion
- Easy-to-read day/night LED display
- Compatible with Solaris line of reusable sensors from Adult to Neonate



Pulse Oximeter Finger Type – iO2

Product Code: PO.F.IO2



- Small and lightweight design
- Dual color OLED display
- Display SpO2, PR, Pulse bar and Plethysmogram
- 3 display modes
- Pulse beat tone function
- Low power consumption, automatic power off
- Low Battery indicator
- Takes 2 AAA alkaline or rechargeable batteries
- Accommodates wide range of finger sizes from pediatric to adult
- CE certified

Pulse Oximeter Finger Type – 500A

Product Code: PO.F.500A

- Small and lightweight design
- Ideal for storing in your jump bag
- Dual colour display
- SPO2 and PR display, PR waveform and bar graph display
- 4 directional display with 8 display modes
- Low power consumption with automatic power off
- Takes 2 AAA batteries and can run for 30 hours
- Alarm function



VS2000 Desktop Vital Signs Monitor

Product Code: VSM.VS2000.CONFIG.A



- The VS200 desktop Vital Signs Monitor is an excellent choice for both EMS and the Hospital Ward. With a range of parameters to choice from this portable unit can perform, NIBP, Pulse Oximetry, EtCO2 (optional), Temp (optional) and 5 lead ECG monitoring (optional)
- Easy to read 7 inch LCD Display
- Weighs only 2kg
- Accurate SPO2 measurement in low perfusion states
- Multiple monitoring parameters can be added to unit
- Auto power off function for saving battery life
- Audible and visible alarms
- Long battery life
- Suitable for adult, paediatric and neonatal patients
- Excellent range of accessories including full range of SpO2 probes, airway adapters, BP cuffs and mounting options (including portable stand)
- Reliable performance
- CE and ISO 13485 Certified



Volumetric Infusion Pump – AJ5808

Product Code: VOL.IN.PMP.5808

- Aluminium alloy casing
- Durable and rust proof
- Linear peristaltic pumping mechanism
- Compatible with most IV sets
- Infused volume range 1-9999ml
- Flow-rate range 1-1200ml/hr. in 1ml/hr. increments
- 3 hours of continuous operation at 25 ml/hr. with built-in rechargeable Ni-MH battery



Fluid Warmer – FW300

Product Code: VOL.IN.W.FW300



- Warms up fluids via an integrated heating plate
- Fluid warmed to between 35 – 40°C
- Warms up IV fluids and enteral nutrition
- Relieves pain associated with cold fluids
- Reduces the potential for hypothermia
- Can be powered directly by AJ5808 infusion pump or use separately with adapter
- Compatible with all administration sets
- Overheating protection

Ambulatory Syringe Pump – CA-500

Product Code: SYR.DVR.CA500

- Accurate, reliable infusion solution for Thalassemia treatment, pain management, chemotherapy, epidural & neonatal care
- Bright OLED Screen
- Treatments with low rate small volume
- Compatible with any brand syringe from 1 to 20ml
- Programmable in ml/hr.
- Powered by 3 AA batteries
- Carrying Bag



Syringe Pump – AJ5803

Product Code: SYR.DVR.ANG03



- Micro infusion syringe pump
- Wide flow rate range from 0.1 ml/h to 1500 ml/h
- Compatible with syringe sized from 10 - 100 ml
- Automatic accuracy calibration
- Adjustable occlusion pressure
- Operates from AC and internal battery (3 hours run-time @25ml/hr.)
- Net mass: 1.85kg
- 2 year warrantee

Portable Syringe Pump– AJ5805

Product Code: SYR.DVR.ANG05

- Lightweight and compact
- Bright OLED Screen
- Compatible with 1 – 20ml syringes
- Infusion rate of 1 – 99ml/hr
- Programmable bolus function
- PCA function – patient can self-administer additional bolus for pain relief
- Carry bag



Portable Syringe Pump – AJ5807

Product Code: SYR.DVR.ANG07



- AC 220V or 110V; built-in rechargeable Ni-MH batteries for emergency
- Power consumption less than 6VA
- 20 ml, 50/60 ml, 100 ml disposable syringes; Compatible with major brand syringes
- Time Settings - 4min, 8min, 12min, 16min and 20min
- Infusion Accuracy ± 20 sec
- Flow rate - 60 ml/h - 1500 ml/h
- Automatic calculation of the delivery rate based on and time pre-selection

Diagnostic ENT Set

Product Code: DIAG.SET

- Chrome-plated bent arm illuminator for lighting the oral cavity
- Two chrome-plated laryngeal mirrors
- Ophthalmoscope with lens wheel and correcting lens diopters ranging from -20 to +20
- Chrome-plated nasal speculum
- Otoscope supplied with 3 reusable specula
- Swivel magnifying lens
- Tongue blade holder
- Handle operated from 2 x C batteries (not included)
- Packaged in reusable high quality box



Additional ENT Diagnostic sets will be added shortly

Single Head Stethoscope

Product Code: GLUC.PLUS



- Economical model
- Aluminium one-sided chest piece
- Screw on metal retaining ring
- Chrome-plated external spring binaurals
- Quality single lumen tubing
- Latex free
- Available in various colours
- Prod Ref: SF100b

Single Head Stethoscope - Deluxe

Product Code: STETH.SH.DELUX

- Features a uniquely designed chest piece for both low and high frequency auscultation
- The Nickel-plated chest piece is finger-contoured for an easy-grip and secure hold
- Large sized Internal spring binaurals
- Non-chill rim design provides patient comfort
- Thick single-lumen tubing
- Accessory box includes one pair of firm ear tips ,one pair of soft ear tips and one diaphragm
- Latex free
- Prod Ref: SF402



Dual Head Stethoscope - Adult

Product Code: STETH.DUAL.DR



- Features an Aluminum dual- sided chest piece
- Non-chill ring
- Screw-on metal retaining ring
- Chrome-plated external spring binaurals
- Quality single lumen tubing
- Latex free
- Available in various colours
- Prod Ref: SF200

Dual Head Stethoscope - Adult

Product Code: STETH.DUAL.DR

- Features a smaller sized chest piece for paediatric auscultation
- Quality single lumen tubing
- Chrome-plated external spring binaurals
- Non-chill ring
- Latex free
- Available in various colours
- Prod Ref: SF201



Dual Head Stethoscope - Lightweight

Product Code: STETH.DH.LW



- Features a lightweight aluminium dual-sided chest piece
- Non-chill retaining ring and bell ring
- Chrome-plated internal spring binaurals
- Thick single lumen tubing
Latex free
- Available in black only
- Prod Ref: SF202

Foetal Stethoscope

Product Code: STETH.FOETAL

- Designed for superior auscultation of the fetal heartbeat
- Chrome-plated internal spring binaurals
- Thick single lumen tubing
- Vinyl coated headrest
Latex free
- Prod Ref: SF207



Stethoscope - Sprague Rappaport

Product Code: STETH.RAPP

- The most popular stethoscope on the market today
- Features 5 interchangeable chest piece fittings to offer the most complete diagnostic capability
- 22" thick double tubing minimises extraneous noise
- Available in black, blue, burgundy, green and pink
- Accessory box includes one pair of mushroom ear tips, one pair of firm ear tips, two diaphragms and three sizes of bells
- Latex free
- Prod Ref: SF301



Dual Head Stethoscope - Deluxe

Product Code: STETH.DH.DLX



- Features a new design chrome-plated chest piece
- Non-chill retaining ring and bell ring
- Quality single lumen tubing
- Internal spring binaurals with plastic soft ear tips
- Latex free
- Prod Ref: SF4011

Stethoscope - Classic

Product Code: STETH.CLASSIC

- Features a stainless steel dual head chest piece for superior acoustical performance
- Stainless steel binaurals
- Non-chill rings for patient comfort
- Quality plastic ear tips ensure a comfortable fit
- Single lumen PVC tubing with thick wall
- Accessory box includes one pair of firm ear tips ,one pair of soft ear tips and one diaphragm
- Latex free
- Prod Ref: SF502



Stethoscope - Classic

Product Code: STETH.CLASSIC



- Approximately 28" overall length
- High acoustic sensitivity for quality auscultation
- Solid Stainless steel chest piece and binaurals
- 2-sided bell & diaphragm combo chest piece
- Cardiology style tubing with internal binaural spring
- Non-chill rim design provides patient comfort
- Accessory box includes one pair of firm ear tips ,one pair of soft ear tips and one diaphragm
- Latex free
- Prod Ref: SF501

Drug Test – Multi 5 Panel

Product Code: DRUG.MUL.5P

- Requires urine sample only
- Excellent stability and accuracy
- Tests for:
 - Amphetamines
 - Methamphetamines
 - Morphine
 - Cocaine
 - Tetrahydrocannabinol (Marijuana)



Drug Test – Multi 5 Panel

Product Code: DRUG.MUL.5P



- Requires urine sample
- Tests for:
 - Opiates
 - Phencyclidine (PCP)
 - Methamphetamines
 - Cocaine
 - Tetrahydrocannabinol (Marijuana)
 - Morphine

Glucometer –Glucoplus

Product Code: GLUC.PLUS

- Requires small blood sample, just 1.5µl
- Accurate results in less than 15 seconds
- 192 test memory with date and time
- Large display



[Back to Contents Page](#)

Disinfectants

Be Safe Hand Sanitiser

Product Code: HS.30ML

- Ingredients: 83% ethyl alcohol; 0.5% chlorhexidine gluconate and selected emollients(Suitable for hands)
- Great as a hand disinfectant to reduce cross infection
- Available in 30, 125, and 500ml and 1, 2, and 5L sizes



Distel High Level Medical Surface Disinfectant

Product Code: D.DIST.NF



- Distel formulation is based on nano-particulate technology enabling the active ingredients to be carried rapidly through cell walls of micro-organisms. The synergistic blend of active ingredients produces the fastest acting and most efficient virucidal and sporicidal disinfectant available in non-oxidising technology
- Distel High Level Medical Surface Disinfectant is suitable for all surfaces, furniture and appliances within hospital, EMS, dental and general practice environments
- Distel High Level Medical Surface Disinfectant can be diluted to various concentrations depending on application
- Available in various volumes
- Available in Eucalyptus and in non-scented
- Also available in wipes (70's)

Medistel Instrument Disinfectant

Product Code: D.MEDI.INSTR.DIS.1LT

- Medistel Instrument and Equipment Disinfectant dissolves and remove proteinaceous and organic material
- Combines secure microbiological credentials with safety in use and a quick contact time
- Medistel is suitable for both human and veterinary applications for high-level disinfection of surgical instruments and medical devices at a dilution of 1:20
- Medistel disinfects in 10 minutes and sterilises in 30minutes
- Available in 1L bottles



OdoGene (Odostel) Odour Eliminator

Product Code: D.ODO.200ML



- OdoGene Odour (soon to be renamed Odostel) Eliminator/Re-Odouriser is based on a powerful blend of essential oils and bactericides
- It is water-based and safe to use on all surfaces except polished wood
- Available in 200 and 500ml and 1, 5 and 20L sizes

Perascope (Perastel) Multishot Product Code: D.PER.7.5LT

- PeraScope (soon to be renamed Perastel) is a new instrument steriliser for heat labile endoscopes, which combines secure microbiological credentials with safety in use and a rapid activity
- Tests conducted in house and with leading manufacturers have shown that Perascope will not harm sensitive instruments including flexible fibre optic instruments
- Effective in 2 minutes against spores, mycobacteria, viruses, fungi and bacteria



[Back to Contents Page](#)

Disposables

Absorbent Granules - Body Fluid Spills

Product Code: ABS.G.100G

- These absorbent granules assist in the speedy removal of body fluid spills and reduce the risk of cross contaminating the surrounding areas
- They are non-toxic, non-corrosive and non-bleaching and can be used on practically any surface
- Suitable for cleaning up blood, urine and emesis spills
- Once the granules are placed on a spill, they react with the spill contents to form a gel. The gel can be easily scooped up and disposed of thereby eliminating the spillage
- Available in 1, 50 and 100g bottles



OdoGene (Odostel) Odour Eliminator

Product Code: D.ODO.200ML



- Used to disinfect skin prior to injection or venipuncture
- Non-woven pad
- Saturated with 0.5ml of 70% Isopropyl alcohol
- Individually sealed in laminated foil packaging
- Available in singles or boxes of 200 swabs

Bag - Cotton Ice– 22.5cm Product Code: BAG.CI

- 22.5cm cotton bag with clip on lid
- Reusable
- Place ice blocks or crushed ice in the bag and apply to affected anatomy to reduce and minimise swelling in sports injuries



Bandage - Conforming Product Code: B.100



- High stretch, lightweight bandage for holding dressings in place
- Cotton-nylon, fray resistant construction
- Available in sizes: 25mm X 2m, 50mm X 4.5m, 75mm X 4.5m, 100mm X 4.5m and 150mm X 4.5m
- Individually wrapped

Bandage - Cohesive

Product Code: B.CPC

- Available in sizes: 25mm, 50mm, 75mm and 100mm
- Dry to the touch
- Co-adhesive high bond formula
- 100% breathable with high tensile strength
- Excellent compression properties.
- Non adhesive, yet sticks to itself



Bandage – Elastic Adhesive Sports Wrap

Product Code: B.EAB.EC



- Manufactured from quality materials
- Skin friendly adhesive
- Excellent for supporting strains and sprains
- Available in sizes: 50, 75 and 100mm X 4.5m
- Available in various brands

Bandage - First Aid Dressings

Product Code: FAD

- Range of First Aid Dressings manufactured from high quality cotton
- Available in:
 - FAD 1 – 25x50x1m
 - FAD 2 – 50x75x1m
 - FAD 3 – 75x100x2.25m
 - FAD 4 – 150x200x4m
 - FAD 5 – 150x200x2.25m
 - FAD 6 – 200x300x6.5m
- Supplied with wrap packaging and number



Kinesio Sports Tape

Product Code: TAPE.KIN.7X4.W



- A lightweight flexible tape that provides versatile support and that promotes better circulation, better support and prevention of joint and muscle injuries
- Suitable for use on ankles, knees, thighs, forearm and wrist and shoulders and upper back
- Caution: wrapping too tightly will impede circulation and is not desirable

CritiBand Multi-Trauma Compression Bandage

Product Code: B.CRIT.100.C

- CritiBand is the latest multi-trauma bandage developed and manufactured by Be Safe Paramedical
- Specifically designed for those instances where maximum direct wound pressure is required to control severe external haemorrhaging
- Designed with the end-user in mind, i.e. can easily be self-applied
- Stitched in trauma pad
- Super stretch - for controlled application of direct wound pressure
- Self-securing - no need for pins or tape, CritiBand wraps onto itself
- Due to its generous length, CritiBand can be applied equally effectively to the head, torso, abdomen and extremities
- Due to its broad structure, pressure is applied in a diffuse manner and therefore does not create a tourniquet effect



Burn-Eaz® Dressings

Product Code: BE.10X10S



- Developed and manufactured by Be Safe Paramedical
- A synergistic combination of Aloe Ferox and Tee Tree Oil brings a new dimension to the treatment of burn injuries
- Excellent antiseptic and heat dissipating properties
- Provides an antiseptic covering and promotes tissue healing
- Individually wrapped in a foil pouch
- Burn-Eaz® dressing are available in sizes: 10 X 10cm, 20 X 20cm and 60 X 40cm
- Sterile for single patient use
- Also available in gel sachet with/without dressing

Burn-Eaz® Gel

Product Code: BE.PP.GEL

- 100g of Burn-Eaz® Gel in a tube
- Developed and manufactured by Be Safe Paramedical
- A synergistic combination of Aloe Ferox and Tee Tree Oil brings a new dimension to the treatment of burn injuries
- Excellent antiseptic and heat dissipating properties
- Apply liberal amounts of Burn-Eaz™ Gel to affected area for instant relief
- Great for soothing the discomfort associated with sunburn



Burn-Eaz® Gel

Product Code: BE.PP.GEL



- 125ml of Burn-Eaz® solution in a handy spray bottle
- Developed and manufactured by Be Safe Paramedical
- A synergistic combination of Aloe Ferox and Tee Tree Oil brings a new dimension to the treatment of burn injuries
- Excellent antiseptic and heat dissipating properties
- Spray liberal amounts of Burn-Eaz™ solution onto affected area. Reapply as needed
- Excellent for sun-burn
- Avoid spraying into eyes

Butterfly Scalp Vein Needle

Product Code: NDLS.BUT

- High quality stainless steel bevelled needle
- Needle protector manufactured from medical quality PE
- Double wings manufactured from medical grade PVC
- 30cm long non-pyrogenic, non-toxic PVC tubing
- Available in 18G, 19G, 21G, 22G, 23G and 25G



Cotton Wool

Product Code: CW.500G



- Highly absorbent
- 100% pure cotton wool
- Suitable for all wound dressing applications
- Excellent for support undercast for bandages
- Available in 50, 100, 250 and 500g packaging

Cotton Wool

Product Code: CW.500G

- Highly absorbent
- 100% pure cotton wool
- Suitable for wound cleansing and for cosmetic use
- Available in 5 and 500g packs



CritiPack™: Intraosseous Vascular Access Pack

Product Code: CRITI.P.EI.I



- The CritiPack™ Intraosseous Vascular Access Pack provides all the kit required to perform an intraosseous infusion in the field
- Sterile Nitrile gloves
- 18G intraosseous needle
- 3-way stop cock
- Antiseptic cleanser
- 20ml syringe
- NaCl flush
- Sterile drape
- Securing Tape
- Tick box skills sheet included

CritiPack™: Maternity and New Born Care Pack

Product Code: CRITI.P.M.I

- The CritiPack™ Maternity and New born Care Pack provides all the kit required to perform a successful emergency delivery in the field
- Sterile drape
- Sterile Nitrile gloves
- Cord clamps
- Disposable scalpel
- Bulb syringe
- Maternity pads
- Gauze pads
- ID bracelets for mom and new born
- Tick box skills sheet included



CritiPack™: Nasogastric Tube Pack

Product Code: CRITI.NG.T.I



- The CritiPack™ Nasogastric Tube Pack provides all the kit required to perform nasogastric or orogastric tube placement in the field
- Sterile Nitrile gloves
- Sterile lubricant
- Protection sheet
- 14 & 16FG nasogastric tubes
- Collection bowl
- Collection bag
- pH paper
- Securing tape
- Tick box skills sheet included

CritiPack™: Needle Crico & Thoracostomy Pack

Product Code: CRITI.P.NC.I

- The CritiPack™ Needle Cricothyroidotomy/ Thoracostomy Pack provides all the kit required to perform needle cricothyroidotomy and thoracostomy in the field
- Sterile Nitrile gloves
- Disinfectant pads
- Oxygen tubing
- 2 x 14G over-the-needle IV cannulas
- 20ml syringe
- Be Safe Paramedical CritiJet Jet Insufflation Device
- Be Safe Paramedical Critivalve (one-way valve)
- Securing tape
- Tick box skills sheet included



CritiPack™: Surgical Airway Pack

Product Code: CRITI.P.SA.I



- The CritiPack™ Surgical Airway Pack provides all the kit required to perform a surgical airway in the field
- Sterile Nitrile gloves
- Disinfectant pads
- 5.5 & 6.0 cuffed ET Tubes
- 20ml syringe
- Catheter mount
- Disposable scalpel
- Sterile pads
- Securing tape
- Tick box skills sheet included

CritiPack™: Suture Pack

Product Code: CRITI.P.MST.I

- The CritiPack™ Suture Pack provides all the kit required to suture uncomplicated lacerations in the field
- Protection sheet
- Surgical gloves
- Plastic bowl
- Cotton wool balls
- Folded Kleenex medical towel
- Gauze swabs
- Antiseptic solution
- 3.0 Blue Monofilament Suture
- 5.0 Blue Monofilament Suture
- Stainless steel needle holder
- Disposable tissue forceps
- Tick box skills sheet included



CritiPack™: Umbilical Vein Cannulation Pack

Product Code: CRITI.P.UV.I



- The CritiPack™ Suture Pack provides all the kit required umbilical vein cannulation on a neonate in the field
- Protection sheet
- Sterile Nitrile gloves
- Disinfectant pads
- Disposable scalpel
- Antiseptic solution
- Umbilical vein catheter size 5
- 20ml syringe
- Saline flush
- Securing tape
- Tick box skills sheet included

CritiPack™: Urinary Catherisation Pack

Product Code: CRITI.P.UC.I

- The CritiPack™ Urinary Catheter Pack provides all equipment required to perform urinary catheterisation in the field
- Protection sheet
- Surgical gloves
- Cotton wool balls
- Antiseptic solution
- Urinary catheters size 16 & 20FG
- Urine bag with tubing
- Sterile lubricant
- Sterile water
- 20ml syringe
- Tick box skills sheet included



CritiSeal – Chest Wound Seal

Product Code: CRITI.SEAL



- Developed and manufactured by Be Safe Paramedical for the treatment of Open Pneumothorax
- Excellent pad adhesion to chest irrespective of presence of hair moisture and blood
- One-way valve promotes thoracic cavity venting on expiration
- One-way valve stops air from being drawn into thoracic cavity during inspiration
- Large valve surface area reduces the chance of clogging due to blood coagulation
- Thick gel adhesive supports an excellent seal and superior chest wall adhesion
- Packaged individually in sealed foil pouch – waterproof and with excellent heat reflecting properties
- Pack includes a 100x100 wound dressing
- The most competitively priced chest wound seal on the market



Electrolyte Rehydration Powder

Product Code: ELEC.P

- Electrolyte replenishment/ supplement
- 5.48g of powder for mixing a 200ml solution
- The formula of Electropak is as recommended in the S.A.M.J.1.
- When correctly reconstituted Electropak solution is isotonic
- Electropak is beneficial for anybody who is dehydrated e.g. diarrhoea, gastro-enteritis and sports people (but not for people who are dehydrated enough to require intravenous therapy)
- The glucose in Electropak is rapidly absorbed, creating the necessary osmotic pressure differential between the intestine and plasma to optimise the water and electrolyte uptake



Eye Shield - Pink

Product Code: EYE.SH



- Manufactured from quality plastic
- Head band for comfy fit
- Available in singles

Sterile Eye Pad

Product Code: EYE.PAD.STER

- Single use sterile eye pad
- Manufactured from soft cotton
- Individually packed



Sterile Eye Pad

Product Code: EYE.PAD.STER



- Sterile eye pad sewn into a conforming bandage
- Pad cushions the eye and absorbs any exudate
- High quality stretch bandage

Gauze Swabs

Product Code: GAU.SWABS

- Available in sterile and non-sterile options
- Excellent for packing wounds
- Soft and durable with a high levels of absorbency
- Low linting
- Produced and packed in clean conditions to ensure total reliability
- Available in various sizes



Hypodermic Needles

Product Code: NDLS.



- Sterile and individually packed
- Manufactured from medical grade stainless steel
- Superbly sharp point minimises pain during puncture
- Colour coded for easy identification
- Available in size 30 to 15G

Identity Bands

Product Code: ID.BAND

- Manufactured from soft PVC with clip-over system
- Pre-printed fields: name, bed, sex and date
- Available in adult and paediatric sizes
- Available in colours: white, yellow, pink, green, blue and red



Intraosseus Needle

Product Code: INTRAOS.NDL



- Manufactured from medical stainless steel
- Superiorly sharp point
- Adjustable depth guard provides depth control for greater patient safety and technique
- Twist-off cap holds the stylet securely in place
- Available in 15, 16 and 18G
- Individually packed
- Single use only

Intravenous (Over-the-needle) Cannula


Product Code: IVC.14G

- Ultra sharp triple edge bevelled stylet
- Plastic catheter
- Sterile
- Individually sealed in tamper evident packaging
- Colour coded according to international standards
- Available in sizes 14 – 24G



Linen Saver

Product Code: LIN.SAV



- Strong and super absorbent Linen saver
- Ideal for overnight use
- 4-ply construction keeps bedding dry & reduces skin irritation
- Excellent value for money

Lubricating Gel

Product Code: LUB.GEL

- Non petroleum based lubricant
- One of the most popular medical lubricants on the market
- Sterile
- Available in 2.5g single sachets and 50 and 100g tubes
- Excellent for lubricating nasopharyngeal airways, endotracheal tubes and urinary catheters



Mask – with loops

Product Code: MASK.SUR.W.LPS



- 3-ply non-woven surgical mask
- Elastic ear loops for comfortable fit
- Breathable for extended application
- Excellent quality

Mask – with tie backs

Product Code: MASK.SUR.TB

- 3-ply non-woven surgical mask
- Adequate tie back length allows for comfortable adjustable fit
- Breathable for extended application
- Excellent quality



Med-Ice Instant Ice Pack

Product Code: MI



- Instant disposable ice pack
- Activates by rupturing internal bag
- Can be kept in first aid / Jump bag and activated when required
- Excellent for applying to sprains, strains, insect bites and stings for instant relief
- Can be applied under arms and in groin area in cases of heat stroke/ exhaustion
- Cover with bandage or cloth to retain cold for extended periods of time

Mop Cap

Product Code: M.C

- Lightweight, non-woven, polypropylene caps
- Lightweight headband secures all hair-styles
- Material is breathable and low linting
- Chemically unreactive and biologically inert



Paper Tape

Product Code: PT.12.5X3 / PT.25X3



- Chemically unreactive and biologically inert
- Non-woven Paper
- Latex-Free
- Hypoallergenic
- Excellent air permeability
- Extremely easy tear by hand
- Available in various sizes

Plasters

Product Code: PL.20.F.ASS

- High quality adhesive plasters
- Individually wrapped
- Available in washable and fabric type
- Available in assortment box and singular design options



Rescue Blanket (Space)

Product Code: RB



- Manufactured from aluminised PET film
- Measures 2m x 1.1m
- Reflects 80% of radiated body heat
- When fully opened it is detectable by electronic devices and highly visible from the air

Reusable Handy Warmer

Product Code: REUSE.HW

- Activate heat by compressing internal disk
- To reheat, simply boil in water for 10 minutes until solution is clear and cool. Reactivate by again compressing internal disk
- Excellent for sports injuries, menstrual pain, tooth ache and arthritis pain



Reusable Hot & Cold Pack

Product Code: REUSE.HC.P



- Simply place pack in freezer for 2hrs to use as cold pack
- Does not freeze
- Place pack in hot water for 7 minutes to use as heat pack
- Non-toxic

Spiroderm Transparent Wound dressing

Product Code: SPIRO

- Transparent polyurethane film dressing
- Hypoallergenic acrylate adhesive
- Excellent for securing IV cannulas
- Available in small and large
- Oxygen and moisture permeable



Syringes

Product Code: SYR



- High quality single use sterile syringe
- Manufactured from medical grade PE
- Available in sizes 1, 2.5, 3, 5, 10, 20 and 50ml
- Available with Luer Lock or Luer Slip
- Available with or without needles

Transparent Tape

Product Code: TAPE.ST.12MM

- Manufactured from PE
- Breathable
- Hypoallergenic and gentle on the skin
- Bi-directional tearing
- Available in size: 12, 25, 50 and 75mm



Tongue Depressors

Product Code: TONG.DEP.SNG



- Manufactured from flexible smooth finished wood
- Rounded ends
- Available as singles or in boxes of 100

Tongue Depressors

Product Code: TONG.DEP.SNG

- Manufactured from non-stretchable cotton
- Also available in disposable non-woven
- Non-sterile
- Use on head, appendage or as a secure arm sling



[Back to Contents Page](#)

EMS & Safety Equipment

Ambulance Helmet with Visor

Product Code: AMBU.H.F

- Manufactured from high quality fibreglass
- Excellent shock absorbing qualities
- Padded inner promotes comfort and safety even when worn for long periods of time
- Fully adjustable head strap allows the user complete control over comfort and stability
- Visor attached to helmet and is adjustable (can be slid into helmet when not in use)
- Shape promotes maximum coverage of head and therefore maximum protection
- Neck protection sheath
- Available in a variety of colours



Ambulance Helmet without Visor

Product Code: HELMET.FIBRE



- Manufactured from high quality fibreglass
- Excellent shock absorbing qualities
Padded inner promotes comfort and safety even when worn for extended periods of time
- Fully adjustable head strap allows the user complete control over comfort and
- This option is suitable for basic head protection on emergency scenes
- Available in a variety of colours

Ambulance Helmet without Visor

Product Code: HELMET.FIBRE

- Range of quality badges representing various EMS qualifications
- Star of Life badge available in Blue and Red
- Flight Medic badge available in blue and red
- Name badge available in plastic or metal, with pin or magnetic backing
- Vinyl reflective Star of Life also available in various sizes



Belt Cutter/ Window Punch

Product Code: BC.WP



- Window punch constructed from solid stainless steel
- Handle manufactured from hard plastic
- Seat belt cutter blade manufactured from stainless steel
- Seat belt cutter blade located in groove in handle. The blade is therefore not exposed and poses no threat to the user
- Channels on either side of blade promote quick application
- Luminous light activated triangle on handle which produces a bright green glow which aids in locating tool should it be dropped in dark place

Emergency Rescue Knife

Product Code: KNIFE.RES.SER

- This Rescue Knife features a 8.6cm drop point blade made of 440 Stainless Steel
- The blade has a non-reflective matt finish and dual thumb studs for ambidextrous one hand opening
- The contoured anodized aluminium handle measures 12.1cm in the closed position and has a blue and white Star of Life logo shield inlay
- The knife sports a seatbelt cutter, a glass breaker and a liner lock
- A metal pocket clip is mounted at the pivot end of the knife
- Available in smooth or serrated blade



Eye – Protective Glasses

Product Code: EYE.PG



- Glasses manufactured from high impact resistant plastic
- Broad nose band promotes comfort and stability
- Side panels ensure that no secretions/ residue can enter from the side (when head is turned)
- Extended cover (superior aspect of glasses) stops secretions/ residue from entering eyes from above
- Arms are fully adjustable and ensure that the glasses fit comfortably over ears

Eye Protective Glasses – Interchangeable Lenses

Product Code: EYE.PG.ICL

- High grade PC frame
- High grade impact resistant Bayer polycarbonate lens
- EN166 And ANSI Z87.1 (+) certified
- ISO 9001 certified
- Interchangeable lens – clear and brown (99.9% UV block) supplied as set
- Comfortable rubber non-slip nose pad
- Modern and comfortable design
- Supplied with hard carry case



Knee Pads – Hard Foam

Product Code: KNEE.PDS.HRD.F



- Knee pads consist of hard outer plastic shell, backed by foam
- The streamlined design makes it comfortable to wear either outside or underneath work clothes
- Knee pads are designed to fit bigger knees too
- Knee pad straps are fully adjustable and contain enough length to be secured to bulkier frames
- Ridges on external plastic frame promote non-slippage

M.V.A. Jacket

Product Code: MVA.ADJUST

- Manufactured from high visibility, high quality nylon
- Available in a selection of colours and designs
- Reflective star of life located on right front
- 2 large pocket with zips and two small pockets with Velcro snaps located on front
- large heavy duty front zip
Large star of life can be placed on rear of jacket further improving visibility
- Available in small, medium and large



Ring Cutter

Product Code: 27-220.16



- High quality stainless steel design
- Super abrasive disk which cuts through rings (except Titanium) quickly and safely

Rescue/ Extrication Gloves

Product Code: GLOV.EXT

- Abrasion resistant leather palm
- Towel lined
- Velcro wrist snuggers
- 3M Scotchlite and nylon grip
- Available in S, M and L



Road Cone – Collapsible

Product Code: RC.C.WL



- Collapsible high visibility road cone
- Reflective band around circumference
- Built in light, which is activated by pressing button at base of cone
- Light has 2 settings: flashing and continuous beam
- Compact design makes storage a breeze

Torch – Rechargeable LED

Product Code: TOR.RE.LED

- Manufactured from Impact-resistant rigid aluminium alloy
- LED Rechargeable Torch
- 220V/12V rechargeable Li-ion battery
- Adjustable settings: high beam, low beam, SOS flasher
- LED Life: 100 000 hours
- Battery: 5 – 6 hours continuous use (dependent on beam selected)
- Dimensions (LXWXH): 35 X30 X155mm
- Supplied with AC/DC charger/ Li-ion battery



Torch – Silver Penlight with gauge/ ruler

Product Code: TOR.PEN.SIL



- Quality pupil penlight torch in silver
- Printed pupil gauge
- Printed ruler
- Smooth surface makes cleaning easy
- Excellent lighting and excellent value

Torch – Penlight – Semi-disposable

Product Code: TOR.PEN

- Quality pupil torch manufactured from hard plastic
- Takes 2 AAA batteries
- Printed pupil gauge on body of torch is useful when measuring pupil size in compromised patients
- Provides excellent lighting
- Cost effective enough to dispose of after several patient contacts



[Back to Contents Page](#)

First Aid Equipment

Collar and Cuff

Product Code: COLL.CUFF

- Soft foam with cotton stockinet overlay
- 4m length
- Plastic ties included
- Excellent as a balanced arm sling, wrist support sling and shoulder immobiliser



Hip Holster – Belt Loop (empty)

Product Code: HIP.HOL.BELT.LP



- Manufactured from high quality nylon webbing
- Excellent for storing pupil torch, scissors etc...
- Attaches to belt with Velcro

Hip Holster – Clip on (empty)

Product Code: HH.UCO

- Manufactured from high quality nylon webbing
- Excellent for storing pupil torch, scissors etc...
- Attaches to belt with loop and Velcro
- Equipment not included



Kidney Dish - Plastic

Product Code: KID.D.PL



- Manufactured from high quality plastic
- Autoclavable
- Available in 15, 200 and 300ml volumes

Kidney Dish – Stainless Steel

Product Code: KID.D.SS

- High quality stainless steel kidney bowl
- Smooth surfaced for easy disinfecting
- Available in 15, 20 and 25cm sizes



Pedometer

Product Code: PED.DMC-03



- Highly visible LCD screen
- Step counts from 0 to 99999
- Mileage, kilometer and calorie count
- Automatic power off when not in use
- Low-power consumption

Scissors – 10cm
Product Code: SC.FA.10CM

- 10cm blade
- Solid stainless steel blade
- Sturdy plastic handle
- Convenient angled design
- Excellent for your first aid kit



Scissors - Universal
Product Code: SC.UNIV



- Stainless steel blade
- Sturdy plastic handle
- Blunt tip with snip protection
- Convenient angled design
- Excellent for cutting dressing and bandages
- Available in 16 or 19cm size
- Available in various colours

Thermometer – Digital

Product Code: THERM.DIGITAL

- 1 minute readout
- Dependable accuracy
- Easy to read digital screen
- Stored memory function
- Automatic shut-off
- Available in flexible tip, waterproof and cartoon character for kids option



Thermometer – Ear infrared

Product Code: THERM.EAR.INFRD



- Fast readout
- 1 minute shut-off
- Uses battery power sparingly
- Operates on 3 AAA batteries
- Soft touch keys
- Spare ear probe covers available

Thermometer – Forehead infrared

Product Code: THERM.IR.F

- Digital backlight display
- Range 30 – 50°C
- Response time < 1sec
- Distance < 10cm
- Mass: 180g



Thermometer – Fridge with alarm

Product Code: THERM.FRIDGE



- High temp and low temp settings
- Velcro fastener
- Large digital display
- Accurate and durable

Thermometer – Oral and Rectal

Product Code: THERM.ORAL

- Oral and rectal version available
- Solid glass tube
- Mercury filled channel
- Temp range 35 – 42°C
- Plastic protective case included



Tourniquet – Jet Pull

Product Code: TOUR.JP



- One of the most widely used tourniquets in peripheral vascular access
- Elasticised cotton belt
- ABS plastic buckle
- Quick release mechanism
- Not suitable for controlling arterial haemorrhaging
- Available in black, dark green, navy blue, royal blue, red and orange

Tourniquet – Velcro

Product Code: TOUR.V

- Economical option
- Elasticised cotton
- Quality Velcro loops and hooks
- Excellent for use during venepuncture procedures



Tourniquet – Rotation Compression (RCT)

Product Code: CRIT-TOUR



- One, or two-handed application
- Velcro and buckle securing system
- Built in padding reduces the potential for "buckle pinching"
- Windlass tensioner
- High quality and an extremely competitive price
- CE and ISO13485 quality accredited
- Available in black only

Tourniquet – Rubber Bandage

Product Code: TOUR.R.100



- Manufactured from neoprene
- 100mm x 2.2m
- Excellent traction
- Cost effective and durable

Pelvigrip Pelvic Binder

Product Code: PEL.G.4.C

- Stabilises pelvic fractures
- Manufactured from neoprene which is light, durable and minimises the potential for the development of pressure sores
- Consists of 4 sizes in one kit, S, M, L and XL
- Its design promotes rapid application with minimal pelvic manipulation
- It's 100% X-ray translucent
- Supplied with quality bag with labelled pouches



[Back to Contents Page](#)

First Aid Kits

Boat First Aid Kit – 4-man

Product Code: 110

- Plastic waterproof box or vinyl pouch
- First Aid Scissors (x1)
- Plastic forceps (x1)
- Safety Pins (x3)
- Cotton Wool (x5)
- Large First Aid Dressing (x1)
- Small First Aid Dressing (x1)
- 75mm Bandage (x1)
- 50mm Bandage (x1)
- Disposable Triangular Bandage (x1)
- Plaster Strips (x10)
- CPR Mouthpiece (x1)
- Large Disposable Examination Gloves (x1 pair)
- Headache Tablets (x10)
- Antiseptic Liquid (75ml)



Boat First Aid Kit – 6-man

Product Code: 130

- First Aid Scissors (x1)
- Plastic forceps (x1)
- Safety Pins (x12)
- Cotton Wool (x10)
- Large First Aid Dressing (x5)
- Small First Aid Dressing (x2)
- 75mm Bandage (x2)
- 50mm Bandage (x2)
- Disposable Triangular Bandage (x1)
- Plaster Strips (x20)
- CPR Mouthpiece (x1)
- Large Disposable Examination Gloves (x2 pairs)
- Headache Tablets (x10)
- Antiseptic Liquid (75ml)



Burn-Eaz® Kit – Home

Product Code: BK.H

- 1 Burn-Eaz® Dressing 20 x 20
- 2 Burn- Eaz® Dressing 10 x 10
- 1 Burn- Eaz® gel 100g
- 1 Cotton Bandage 100mm
- 2 Cotton Bandage 50mm
- 1 Scissors 10/13cm
- 1 Safety Pins 12's
- 1 Latex Gloves Large



Burn-Eaz Kit® – Personal

Product Code: BK.P



- 1 Burn- Eaz® Dressing 10 x 10
- 1 Burn- Eaz® Gel 100g
- 1 Cotton Bandage 100mm
- 3 Plaster Strips Fabric
- 1 Scissors 10/13cm
- 1 Safety Pins 3's
- 1 Latex Gloves Large
- 1 Paper Tape 12.5mm x 3m

Burn-Eaz Kit® – Responder

Product Code: BK.R

- 1 Burn- Eaz® Dressing 60 x 40
- 4 Burn- Eaz® Dressing 20 x 20
- 4 Burn- Eaz® Dressing 10 x 10
- 2 Burn- Eaz® Gel 100g
- 1 Burn- Eaz® Spray 125ml
- 6 Cotton Bandages 50mm
- 6 Cotton Bandages 100mm
- 3 Plaster Strips Fabric
- 1 Scissors 16cm
- 1 Pair of Tweezers
- 2 Safety Pins 12's
- 3 Latex Gloves Large
- 3 latex Gloves Medium
- 2 Paper Tape 12.5mm x 3m



Burn-Eaz Kit® – Restaurant

Product Code: BK.LR



- 2 Burn- Eaz® Dressing 20 x 20
- 4 Burn- Eaz® Dressing 10 x 10
- 1 Burn- Eaz® Gel 100g
- 1 Burn- Eaz® Spray 125ml
- 2 Cotton Bandages 50mm
- 2 Cotton Bandages 100mm
- 10 Plaster Strips Blue
- 1 Scissors 10/13cm
- 1 Metal Tweezers
- 2 Safety Pins 12's
- 1 Latex Gloves Large
- 1 latex Gloves Medium
- 1 Paper Tape 12.5mm x 3m
- Kit available in large or small

Emergency Road Side Kit - Large

Product Code: BAG.ER.L

- 1 Warning Triangle
- 1 Set of Jumper Cables
- 1 Tow Rope
- 1 Tyre Weld (Emergency Puncture Repair Kit)
- 1 QE-20 Lubricant Spray (for spark plugs)
- 1 Emergency Light
- 1 Wound Cleaner (75ml)
- 1 Gentian Violet (20ml)
- 1 CPR Mouthpiece
- 1 Pair Disposable Gloves
- 1 Rescue Blanket (foil)
- 1 Set of Metal Forceps
- 1 Pair of 16cm Universal Scissors
- 2 Cotton Wool balls (5g)
- 2 Eye Drop Vials
- 10 Plaster Strips
- 1 Set Safety Pins
- 4 Triangular Bandages
- 1 Burn-Eaz™ Dressing (10x10)
- 1 Gauze Dressings (5)
- 2 First Aid Dressing no. 2
- 1 First Aid Dressing no. 3
- 1 Bandage 50mm



Small kit also available

Extractor Pump Kit

Product Code: EPK



- Kit for extracting venom from stings and bites of all types
- Extractor creates vacuum of 75 – 80 kpa
- Kit includes extractor, various size suctioncaps, plasters, disinfectant swabs and razor

Factory First Aid Kit

Product Code: 109

- 2 x Disposable Gloves Medium
- 2 x Disposable Gloves large
- 100 x Swabs
- 2 x Antiseptic Wound Cleaners
- 8 x Bandages
- 2 x CPR Mouth Pieces
- 4 x Triangular Bandages
- 1 x Elastic Adhesive
- 1 x Plaster Strip
- 12 x Safety Pins
- 2 x Straight Splints
- 10 x Sterile Swabs
- 2 x Cotton wool
- 1 x Metal Forceps
- 1 x Non-Allergenic Adhesive Tape
- 1 x 10cm first aid Scissors
- Available in metal container, plastic box or in nylon bag



Grabber Bag Sports Kit

Product Code: 510GB



- 4 Bandage: Conforming 75mm
- 1 Antiseptic cream - centrimide 25g
- 2 Plaster fabric stretch roll 25 x 3
- 20 Paracetamol Tablets
- 2 Bandage: knit crepe 50mm
- 1 Bandage: conforming 50mm
- 2 Bandage: knit crepe 75mm
- 1 Ung meth sal 25g
- 20 Plaster strips 19 x 72 (fabric)
- 2 Triangular bandage - disposable
- 1 Safety pins 12's
- 1 Skin closure 3 x 75
- 1 Scissors 10/13cm
- 1 Rescue/space blanket
- 1 Paper tape 25 x 3
- 1 Neck collar - grey trauma foam – medium
- 1 Metal tweezer
- 2 Medice - instant ice pack
- 2 Pair - latex gloves large
- 5 Sterile gauze swabs 75 x 75 x 8ply
- 1 Cotton wool 100g
- 1 Burn-Eaz™ dressing 10cm x 10cm
- 2 Bandage - elastic adhesive 75mm
- 2 Bandage - elastic adhesive 50mm

Gun/ Stab Wound Kit

Product Code: 210GS.V

- 2 No. 2 First Aid Dressings
- 2 No. 3 First Aid Dressings
- 2 No. 4 First Aid Dressings
- 1 Chest wound seal
- 1 Scissors 16cm universal
- 1 75ml Savlon wound cleaner
- 1 Cotton wool balls 5g
- 1 Gauze swabs 75x75x8 ply (5's)
- 1 CPR mouthpiece
- 1 Pair of gloves medium
- 1 Pair of gloves large
- 2 Triangular bandages
- 1 Stretch Fabric Plaster 25mmx3m



Be Safe Paramedical IFAK (Individual First Aid Kit)

Product Code: REF:FAK.TACT.P



- Selection of bags to choose from
- Contains:
 - 1 x CritiBand Trauma Bandage
 - 2 x No. 3 first aid dressings
 - 1 x Tape
 - 1 x Rotation Compression Tourniquet
 - 1 x 75mm Gauze Roll
 - 2 x Nitrile Gloves
 - 1 x CPR Mouthpiece
 - 1 x CritiSeal Chest Wound Seal
 - 1 x 28 Fr Nasopharyngeal Airway
 - 1 x Lubricant Sachet
 - 1 x 14G IV needle
 - 1 x Check-Valve
- Custom items as per client's request

Motorist First Aid Kit)

Product Code: 107

- Excellent high quality compact bag which can fit under the seat or in the boot of your car
- 1 Antiseptic solution 100ml
- 2 Swabs
- 2 Cotton Wool 50g
- 4 Triangular bandages
- 4 Bandages 75mm
- 1 Elastic adhesive 25mm x 3m
- 1 Adhesive dressing strips 20's
- 4 First aid dressings no 3
- 1 Paper tape 25mm x 3m
- 2 Straight splints
- 2 disposable gloves large
- 2 Disposable gloves medium
- 1 CPR mouthpiece
- 1 Metal forceps
- 1 First aid scissors
- 12 safety pins



Various packaging options available

Office/ School First Aid Kit

Product Code: 160PLO



- 1 Tri-bandage
- 2 Pain/fever tabs 10'S
- 1 Eyedrops 5 ML
- 1 Wound cleaner 75 ML
- 1 Plaster strips 20's
- 1 Cotton wool 50 G
- 1 Burn-Eaz ® 100 mm X 100 mm
- 1 Safety pins 12's
- 1 Metal tweezers
- 1 16cm Universal stainless steel scissors
- 1 First aid dressing NO.5
- 1 First aid dressing NO.3
- 1 First aid dressing NO.2
- 2 bandages 100 mm
- 2 Bandages 75 mm
- 2 bandages 50 mm
- 1 Gauze 75x75x8 ply (5'S)
- 1 Med Ice
- 1 CPR mouthpiece
- 1 Pair gloves
- 1 Elastic Adhesive 25 mm X 3m

Various packaging options available

Ouch Pouch – Just for you

Product Code: OP.JFY

- 2 saline eyedrops
- 1 medical info sheet
- 1 Antiseptic wipes
- 1 Paracetamol (10's)
- 3 Safety Pins
- 5 Fabric plaster strips 19 x 72 mm
- 1 Disposable pair of gloves
- 3 Cotton buds
- 1 CPR Mouthpiece



Restaurant First Aid Kit

Product Code: 170



- 2 disposable Tri-bandage
- 20 Paracetamol tablets
- 1 Artificial Tears (5's)
- 1 Savlon 75 ML
- 1 Plaster strips (100's)
- 1 Cotton wool 100 G
- 4 Burn-Eaz ® 100 mm X 100 mm
- 2 Burn-Eaz ® 200 mm X 200 mm
- 1 Safety pins 12's
- 1 Metal tweezers
- 1 10cm first aid scissors
- 2 First aid dressing NO.3
- 2 bandages 50 mm
- 1 Gauze 100
- 2 Medlce
- 1 CPR mouthpiece
- 2 Pairs latex gloves (large)
- 2 Pairs Latex Gloves (Medium)

School/ Club First Aid Kit

Product Code: 180

- 2 Tri-bandage
- 1 Wound cleaner 75 ML
- 1 Plaster strips 20's
- 1 Cotton wool 50 G
- 1 Burn-Eaz ® 100 mm X 100 mm
- 1 Safety pins 12's
- 1 Metal tweezers
- 1 10cm first aid scissors
- 2 First aid dressing NO.5
- 2 First aid dressing NO.3
- 2 First aid dressing NO.2
- 2 conforming bandages 100 mm
- 2 conforming bandages 75 mm
- 2 conforming bandages 50 mm
- 1 Elastic adhesive 25 x 3m
- 10 antiseptic wipes
- 1 Plaster Strips 20's
- 2 CPR mouthpiece
- 2 Pairs latex gloves (1 large +1 medium)
- 1 Skin closure
- 1 Rescue Blanket
- 1 Sal Volatile 20ml
- 1 Friars Balsam 20ml



Solas First Aid Kit

Product Code: 100



- 4 Triangular bandages – disposable
- 1 Safety pins 12's
- 1 Scissors - universal 16cm black
- 20 Plaster strips 19x72 – fabric
- 1 Painamol (100's)
- 2 F.A.D no: 6 200x300x6.5m
- 2 F.A.D no: 5 150x200x2.25m
- 2 F.A.D no: 3 75x100x2.25m
- 1 Burn-Eaz ® dressing 200 x 200mm
- 4 Burn Eaz ® dressing 100 x 100mm
- 6 Bandage - conforming 75x4.5
- 2 Antiseptic cream - Cetrimide 50g
- 1 Pair Gloves - latex (large)
- 1 CPR mouthpiece

Sports First Aid Kit Moonbag

Product Code: 520

- 1 Med Ice instant ice pack
- 1 Elastic adhesive 50mm
- 1 Crepe bandage 75mm
- 1 Cotton Wool balls 5g
- 1 Pair latex gloves (Medium)
- 1 Savlon 75ml
- 1 Metal tweezer
- 1 Needle – hypodermic 21G
- 1 pain/ fever tabs 10's
- 1 PVP cream 25g
- 1 Plasters – fabric strips 10's
- 1 Triangular bandage
- 1 Scissors 16cm universal
- 1 CPR mouthpiece
- 2 Eye drops



Immobilisation Devices

Adjustable Stiff Neck Collar

Product Code: NECK.STF.ADJ.A

- 580 x 150mm
- Manufactured from quality HDPE material
- High density polyethylene foam inner
- Hook and loop fastener
- Adjustable to various sizes
- Various colour options available on larger orders
- Also available in paediatric size



Extrication Device - Kendrick

Product Code: EXT.DEV.KED



- Manufactured from heavy duty vinyl coated nylon
- Built in handles
- Colour-coded, sewn-in securing straps and snap lock buckles for quick and easy application
- Includes wrap-around vest, adjustable padded neck roll, two head straps and carrying case
- Multiple size adjustable

Also available in paediatric size

CritiStrap – Reusable Splint Strapping

Product Code: EXT.DEV.KED

- Save costs on traditional flannel strips used for securing Kramer Wire Splints
- CritiStrap is manufactured from high quality stretch nylon
- Velcro hook and loop securing system
- CritiStrap's generous length provides for the size variations commonly encountered by emergency care personnel
- CritiStrap's design allows for straps to be joined to secure bilateral fractures of the upper and lower extremities
- The broad structure of the straps promotes excellent stability and ultimately better patient comfort
- Reuse to save on your equipment budget
- Available in Large (blue label) for adult lower extremities and Medium (red label) for adult upper extremities and paediatric lower extremities
- Large pack consists of 2 Large and 1 Medium CritiStrap in a quality nylon bag
- Medium Pack consists of 3 medium straps in a quality nylon bag



Head Immobiliser (Head Blocks) - Universal

Product Code: HEAD.IMM



- Manufactured from high density foam with plastic coating
- Waterproof and non-absorbable
- Large ear inspection holes
- Can be attached to a scoop stretcher or a spinal board
- Adjustable head straps
- X-ray translucent and CT/ MRI scanning compatible

Head Immobiliser (Head Blocks) – Universal Paediatric

Product Code: HEAD.IMM.CHLD

- Poly dipped vinyl coated foam
- Specifically designed for the paediatric patient, taking into consideration factors such as: size, anatomical differences and facial bulk
- Anatomically correct ear inspection holes
- Base plate can be secured to scoop stretchers and spinal boards
- Dual shaped blocks that can be reversed to fit concave and flat boards
- Adjustable head straps
- X-ray translucent and CT/ MRI scanning compatible



Neck Collar - Rigid

Product Code: NECK.STF



- Manufactured from quality HDPE material
- Flip chin piece
- Available in small, medium, large and extra-large
- Loop contact closures
- Opening on anterior section allows for inspection of trachea and external jugular veins
- Radiotranslucent (X-Ray, CT & MRI compatible)

Neck Collar – Grey Foam

Product Code: NECK.GT

- Manufactured from dense foam which is rigid for stability and yet comfortable to wear for extended periods of time
- Opening on anterior section allows for inspection of trachea
- Loop contact closures
- Available in small, medium and large



Splint - EVA

Product Code: SPL.EVA



- Reusable
- Lightweight and compact
- Radiolucent
- Fastens in place with tape or wrap of choice; no extra equipment is needed
- Not affected by extreme temperatures or altitudes
- Closed-cell foam allows easy cleaning and disinfection; material compatible with all standard cleaning solutions
- Latex Free

Kramer Wire Splint - Padded

Product Code: SPL.K.L

- One of the most widely used splints in the emergency services
- Mesh inner, covered with dense foam and plastic
- Can be easily manipulated to fit angulated fractures
- Plastic outer makes cleaning easy
- Available in 'Long' size for lower extremities and 'short' for upper extremity fractures
- Excellent stability when used in conjunction with Be Safe Paramedical's CriteStraps



Splint – Lock In

Product Code: SPL.PL.GREEN



- Manufactured from solid plastic
- The design allows for the splints to be interlocked to form a larger splint, thus catering for a wide range of needs
- Quick application
- Reusable
- Easy to clean and keep clean
- Small size makes it easy to store in your jump or first aid Bag

Splint – Unilateral Traction

Product Code: SPL.TR.UNI.A.I

- Manufactured from light and durable aluminium
- Adjustable length
- 4 adjustable nylon leg straps
- Adjustable ankle hitch
- Automatic traction mechanism
- Quality bag included
- Also available in paediatric size



Portable Leg Traction Device

Product Code: SPL.TR.PORT



- Compact bag that only measures 28 x 15cm
- Weighs less than 500g
- Fits easily into your jumpbag
- One size adjustable to fit both adult and paediatric patients
- Aluminium alloy traction pole
- Single person application
- Nylon strapping support
- Extremely competitively priced

Spider Harness – Velcro or Buckle Clip

Product Code: STRAP.SP.BUCK

- Used to secure patient to a spinal board, basket stretcher or scoop stretcher
- Manufactured from quality nylon webbing
- Available with Velcro tie or Plastic Clip Buckles
- Available in black only



[Back to Contents Page](#)

Respiratory Care

Bag-Valve-Mask Resuscitator - Adult

Product Code: RES.SIL.A.I

- Reusable silicone resuscitator
- 1500ml volume
- 2500ml PVC/PE reservoir
- "pop-off" valve 60cmH₂O
- Intake valve set
- Reservoir intake valve set
- Autoclavable up to 134°C



Bag-Valve-Mask Resuscitator – Child

Product Code: RES.SIL.CH



- Reusable silicone resuscitator
- 550ml volume
- 2500ml PVC/PE reservoir
- "pop-off" valve 40cmH₂O
- Intake valve set
- Reservoir intake valve set
- Autoclavable up to 134°C

Bag-Valve-Mask Resuscitator - Infant

Product Code: RES.SIL.PAED

- Reusable silicone resuscitator
- 280ml volume
- 600ml PVC/PE reservoir
- "pop-off" valve 40cmH₂O
- Intake valve set
- Reservoir intake valve set
- Autoclavable up to 134°C



Entonox Regulator with Demand Valve

Product Code: ENT.REG



- For use with entonox cylinders
- Brass inserts
- Demand valve activates when negative pressure applied (inhalation)
- Demand valve attaches directly to mask
- Can be attached to any standard entonox face mask
- Crush proof tubing

Entonox Regulator with Demand Valve - Import

Product Code: ENT.REG.I

- For use with entonox cylinders
- Lightweight aluminum demand valve
- High quality and durable regulator with brass inserts
- Crush proof tubing
- Silicone mask



Nasal Cannula

Product Code: NAS.CAN



- Manufactured from non-toxic medical grade PVC
- 2m tubing
- Individually packaged
- Available in adult, child, and neonatal sizes

Nebulizer Mask – Adult and Paediatric

Product Code: O2MNA

- Manufactured from clear, medical-grade PVC
- Adjustable nasal clip
- Oxygen tubing included
- Equipped with a 6ML medication chamber.
- Available in adult and paediatric sizes



Nebulizer - In-line

Product Code: NEB.IL



- Designed to provide B₂ stimulants to intubated patients or patients requiring assisted ventilation with manual resuscitator
- Includes; medication chamber, catheter mount and oxygen tubing
- Catheter mount attaches to regular 15mm ETT adapter
- Excludes manual resuscitator and endotracheal tube

Nebulizer– RJ-202 Compressor

Product Code: NEB.RJ-202

- Extremely durable and reliable
- Built in carry handle
- Easy to use one button operation
- 2 – 10ml medication capacity
- Low noise operation
- Oil-free air piston



Oxygen Mask – Partial -Rebreather

Product Code: O2MPRA



- Made of clear, non-toxic PVC
- Delivers concentrations of 60 – 90% at 6 – 10L/min flow
- Adjustable nose clip
- Available with or without oxygen tubing
- Available in adult and paediatric sizes

Oxygen Mask – Non-Rebreather Product Code: O2MRA

- Made of clear, non-toxic PVC
- Delivers concentrations of 95 – 100% at 10 - 15L/min flow
- Adjustable nose clip
- Available with or without oxygen tubing
- Available in adult and paediatric sizes



Oxygen Mask – 40-60% Venturi Adult and Paediatric Product Code: O2M.40%.WT.A.I



- Made of clear, non-toxic PVC
- Available in 40 or 60% concentrations
- Adjustable nose clip
- Available with or without oxygen tubing

Oxygen Cylinder – 3L

Product Code: O2.424

- Manufactured from quality aluminium with colour code according to SAPS 019 (black cylinder with a white neck)
- Epoxy coated
- Fitted with 2 pin-index valve
- Mass: 3.6kg



Oxygen Cylinder – 10L

Product Code: O2.1400LT



- Manufactured from quality aluminium with colour code according to SAPS 019 (black cylinder with a white neck)
- Epoxy coated
- Available in 2-pin index and bull-nose configurations
- Mass: 8.2kg

Oxygen Regulator – Dial Flow Pin Index

Product Code: O2R.PI

- Aluminium body with brass inlet retainer, relief valve, piston and rotor flow control hosing
- Silicone "O" rings
- Dial Flow: up to 25l/min
- 3000 PSI



Oxygen Regulator - Bullnose

Product Code: O2RB



- Aluminium body with all brass inner core
- 12 flow rates
- 50 PSI Output
- Dial Flow: 0 - 25l/min
- 3000 PSI Gauge

Peak Flow Meter

Product Code: PEAK.FLOW

- Manufactured from recyclable ABS Plastic
- Accuracy: $\pm 10\text{L/min}$ or 10% of the reading
- Altitude effects: Lowers reading by approx. 5% per 1000m
- Measurement range: 50-800L/min BTPS



Ventilator - VENTH100C Transport Ventilator

Product Code: VENT.H-100C



- Mass: 3kg
- Modes: SIMV, SPONT, SIGH and PEEP
- Tidal Volume: 50 – 1500ml
- Breath rate: 5 – 60 B/min and 1 – 12 for SIMV
- %O₂: 45 or 100%
- PEEP: 0 – 10hpa
- Over pressure relief: <60hpa
- Base flow: 0 – 5L/min
- Supplied in bag

Ventilator – MeVenti Mt2 Transport Ventilator

Product Code: VENT.MT2

- The MeVenti Mt2 is a Portable Emergency Ventilator suitable for the ventilation of infants, children and adult patients via endotracheal tube and face-mask
- High end/ low cost
- Ventilation Modes: OFF, AUT+AST, PSV, CPAP
- Minute Volume: 1 – 16L/min
- Tidal Volume: 15 – 3000ml
- Frequency: 5 – 70 breaths/min
- Pressure Limit: 0 – 50cm H₂O (adjustable)
- Integrated PEEP: 0 – 20cm H₂O
- Integrated mixer: 50 or 100% oxygen, or 25% when compressed air is used
- Audio and visual alarms for apnoea, high pressure, low pressure, gas supply failure and power supply failure
- Internal battery: 10 hours operation



Ventilator – Sirio Baby 200 Emergency Transport Ventilator

Product Code: VENT.SB.200



- Designed for ventilation of new-born and premature neonates
- **Ventilatory Modes:**
 - CPAP
 - IMV
 - IPPV
 - PEEP
 - Manual breath
- **Adjustable Parameters:**
 - Inspiratory time
 - Expiratory time
 - Minimum and maximum airway pressures
 - FiO2 concentration from 21 to 100%
 - Respiratory rate
 - PEEP
- Medical Gas: Medical Oxygen or Compressed Air
- **Measured Parameters:**
 - Airway peak pressure
 - Mean airway pressure
 - Respiratory rate
 -
- Various mounting options
- Manufactured in Italy
- 2 Year guarantee



[Back to Contents Page](#)

Stretchers & Lifting Devices

Basket Stretcher – Stainless Steel

Product Code: STR.BASK.SS.FRM

- Manufactured from stainless steel
- Mass: 22kg
- Load bearing: 350kg
- Can be split in two for easy storage and transport
- Can accommodate a spinal board for added patient stability
- Ideal for rotor wing aircraft operations



Basket Stretcher –Plastic

Product Code: STR.BASK.PLS



- Locally manufactured
- Solid 1-piece frame
- High quality fibreglass construction
- PVC railings
- Loading capacity: 350 Kg

Criti Stair Chair Stretcher

Product Code: STR.CRITI.CH

- Manufactured from high strength aluminium alloys
- Utilised for evacuating patients down stairs in tall buildings
- Single person operation
- Special wheel design provides excellent tracking and stability in confined spaces and in narrow stairwells
- Mass: 11 kg
- Load bearing: 159kg
- Can be folded for easy storage



Emergency Trolley YC-3M

Product Code: STR.E.T



- Powder-coated steel frame
- Height adjustable by means of oil-hydraulic foot pump
- Back adjustable by means of pneumatic spring
- Wheel locking mechanism keeps unit stable when moving patients
- Clip in side brackets
- Load bearing: 159kg

Scoop Stretcher – Plastic folding

Product Code: STR.SCOOP

- Manufactured from high density plastic
- Steel frame and attachment points
- Quick release mechanism for easy assembly and disassembly
- Height adjustable
- Mass: 7kg
- Load bearing: 159kg



Scoop Stretcher – Aluminium Folding

Product Code: STR.SCOOP.F



- Manufactured from aluminium alloy
- Quick release mechanism for easy assembly and disassembly
- Mass: 10kg
- Load bearing: 159kg

Stretcher – Aluminium Folding Spine Board

Product Code: STR.AF.SB

- Lightweight aluminium alloy framework
- Restraining straps included
- Unfolded Size (LxWxH): 85x52x5cm
- Folded Size (LxWxH): 94x50x7cm
- Mass: 7 kg
- Load Bearing: 159 kg



Stretcher – Chair Stair

Product Code: STR.CH/ST



- Lightweight aluminium alloy framework
- Two piece chest and thigh restraints
- Telescopic handle
- Can be folded for easy storage
- Mass: 8.5kg
- Load bearing: 159kg

Stretcher – Foldable Chair Stair (Import)

Product Code: STR.ST/CH.YDC1A12

- Manufactured from high-strength aluminium alloy and polyvinyl chloride material
- Fold-away handles conveniently located at rear of stretcher for easy directional movement
- Fold-away large diameter wheels located at bottom rear support bar
- Can be used as basic fold-away stretcher
- Converts into high back stair chair
- 1 x over-lap and 1 x over-chest restraint
- Back support height: 137cm
- Unfolded size: 190 x 58 x 20cm
- Folded size: 105 x 55 x 20cm
- Mass: 14kg
- Maximum load capacity: 160kg



Stretcher – Double Fold Pole

Product Code: STR.DBL.FLD.POL



- Lightweight aluminium alloy framework
- Double fold mechanism for confined space storage
- Mass: 10kg
- Load bearing: 159kg

Stretcher – Single Fold Pole

Product Code: STR.SFP

- Lightweight aluminium alloy framework
- Single fold mechanism for easy storage
- Mass: 6kg
- Load bearing: 200kg



Stretcher – Foldaway with adjustable back

Product Code: STR.FA.AB



- Lightweight aluminium alloy framework
- Adjustable back
- Folds in half for easy storage
- Patient restraining straps included
- Mass: 8kg
- Load bearing: 159kg

Stretcher – Zip

Product Code: STR.PL.ZIP.I

- Manufactured from quality vinyl
- High quality zip mechanism
- 10 integrated carry handles
- Mass: 1kg
- Load bearing: 159kg



Stretcher – Plastic Spine YSC-J02

Product Code: STR.FA.AB




- Manufactured from high quality PE material
- Pre-stressed for additional weight carrying capabilities
- Non-slip mild grain surface
- 20 large gripping handles
- Mass: 7.5kg
- Load bearing: 159kg

Stretcher – Plastic Spinal

Product Code: STR.PS

- Manufactured from high quality PE material
- Angled edges to aid in log rolling patient
- Tapered design follows body contours ensuring easy maneuverability in confined spaces
- Pre-stressed for additional weight carrying capabilities
- Non-slip mild grain surface
- Ten large gripping handles
- Mass: 7kg
- Load bearing: 159kg



Stretcher – Wooden Spinal I

Product Code: STR.TR.B.L

- Available in pine or marine ply
- Compatible with any standard spider harness and head immobilisation device



Stretcher – SLIX 100 Roll Up Rescue Stretcher

Product Code: STR.SLX100



- Designed to be used with a 2 person Tripod for confined space or other forms of rescue
- Flexible and smooth surface allows easy sliding over any surface
- Strong handles allow it to be maneuvered by several rescuers
- Foot loops and webbing straps secure and allow the stretcher to mould around the casualty
- Vertically and horizontally hauling options
- The lifting straps are colour coded for ease of use
- Stretcher mass: 1.4kg
- Load bearing 200kg

Stretcher – Rescue - YSC-R02
Product Code: STR.RES.YSC-R02

- Manufactured from specialised plastic composite
- Light weight and compact
- Suitable for both horizontal and vertical lifting
- Also suitable for confined space rescue and high angle rescue
- Dimensions:
 - Unfolded: 245x93cm
 - Net Weight: 10kg
- Maximum working load: 160kg
- Supplied with bag for easy storage
- CE&ISO9001/13485 certified



Stretcher – Rescue - YSC-R02
Product Code: STR.RES.YSC-R02



- Constructed from spliced, ramin slats (hardwood) with a tough rot proof canvas exterior
- Highly resistant to stains and fluids
- Strong and supportive head restraints offer stability to the head when vertical lifting is necessary
- Patient restraints provide excellent patient stability in both vertical and horizontal position
- Safe working load: 120kg
- Dimensions: 150 X 100cm
- Mass: 10kg

Stretcher – Self-loading

Product Code: STR.SLOADNG.IM.3GW

- Manufactured from high-strength aluminum alloy material
- Height can be adjusted to 12 positions
- Height for loading into ambulance is 62 – 64cm
- Water-proof foam cushioned mattress
- Backrest and leg components are pneumatically adjustable
- Single person operation
- Leg folding levers located at head and foot of stretcher
- Main castor is 150mm in diameter
- Dimensions: 203×56×92cm
- Maximum load: 250kg



Stretcher – Self-loading

Product Code: STR.SLOADNG.IM



- Manufactured from aluminium alloy
- Extremely light, with excellent rigidity and support
- High density foam mattress with PVC covering
- Broad wheel design makes uneven terrain easier to navigate
- Self-loading
- Automatic locking mechanism
- Load bearing 181kg

Stretcher –YSC11 Self-loading

Product Code: STR.SLOADNG.IM.YSC11

- Extremely cost effective without compromising quality
- Manufactured from high-strength aluminium alloy material
- Height can be adjusted to 12 positions
- Thick-walled oval legs with increased strength
- Water-proof foam cushioned mattress
- Backrest adjustable from 0 – 75o
- Single person operation
- Leg folding levers located and head and foot of stretcher
- Main castor is 125mm in diameter
- Dimensions: 200×62×32cm
- Maximum load: 250kg
- Stretcher mass: 40kg
- Accessories available: Oxygen holder and adjustable IV pole



Stretcher – Standard Ambulance

Product Code: STR.AMBU.ST



- Manufactured from aluminium alloy
- Extremely light, with excellent rigidity and support
- High density foam mattress with PVC cover
- Low frame provides excellent stability on a multitude of surfaces
- Includes mounting bracket
- Load-bearing: 159kg
- Available with orange mattress only

Stretcher – Vacuum Mattress

Product Code: VAC.MATT.COMP

- Manufactured from high quality TPU material
- 8 built in handles to ensure safe and stable handling of patient
- 4 built in patient restraining straps
- Smooth surface promotes hygiene and infection control during cleaning
- Vacuum Mattress is X-ray and CT scan friendly
- Double-chamber pump and bag included



Easy-Move Lateral Patient Transfer System

Product Code: STR.EM



- Quick and reliable application used to transfer patients from one bed to another
- Cost effective when compared to the cost of work related injuries
- Excellent patient comfort
- X-Ray Translucent
- Manufactured from high quality nylon (cover) and PP Board (inner board)
- Folds in half for easy storage
- Dimensions: 170 x 50cm
- Mass: 2.3kg
- Load capacity: ≤ 200 kg
- CE marked

[Back to Contents Page](#)

Vascular Access

CritiPack® Emergency Intraosseous Pack

Product Code: CRITI.P.EI.I

- Provides the healthcare professional with all the equipment required to perform intraosseous vascular access
- The pack contains:
 - 1 x Sterile Nitrile Gloves
 - 1 x Protection Sheet
 - 2 x Disinfectant pads
 - 1 x 18G Intraosseous needle
 - 1 x 3-way stopcock
 - 1 x 20cc syringe
 - 20ml Saline flush
 - Gauze Swabs
 - 1 x Securing tape



CritiPack® Umbilical Vein Catheterisation Pack

Product Code: CRITI.P.UV.I



- CritiPack provides the healthcare professional with all the equipment required to perform an umbilical vein catheterisation in a neonate
- The pack contains:
 - 1 x Sterile Nitrile Gloves
 - 1 x Protection Sheet
 - 2 x Disinfectant pads
 - 1 x Disposable scalpel
 - 2 x Sterile size 5 umbilical catheters
 - 1 x 20cc Syringe
 - 2 x 10ml Saline flush
 - 1 x Cord tie
 - 1 x Securing tape

Administration Sets - Intravenous

Product Code: ADM

- Conforms to International Standards
- Flexible Drip Chamber with air vent filter
- Quality low kink tubing with flow controller for accurate & unrestricted flow
- Special beveled extra strong spike.
- Available in 60, 15, 15/20 and 20 drops/ml
- Luer Lock or Luer Slip connector available
- Individually packed
- Sterile, single use only
- Non Toxic and pyrogen free



Administration Set – 10 drop/ml

Product Code: ADM10



- Transparent and flexible drip chamber with built in air vent and blood filter
- High kink resistant tubing
- Efficient flow regulator for accurate adjustment of Blood flow rate
- Luer Lock connector
- Individually packed
- Sterile, single use only
- Non Toxic and pyrogen free

Administration Sets – Controlled Volume 150ml (Buretrol)

Product Code: ADM.CV

- Contains a metered small-volume chamber designed to limit the amount of solution available to the patient
- Suitable for gravity infusions
- Dosage container 150 ml, with fine graduation
- Injection site (latex) for drug addition
- Automatic shut-off swim valve to prevent air embolism.
- Flexible, transparent drip chamber



Intravenous Cannula

Product Code: IVC.BSPM



- Thin-wall construction permits high flow rates
- Radiopaque with 3 x-ray detectable stripes
- Clear and transparent smooth surfaced catheter
- Back-cut special designed needle
- Hydrophobic blood stopper
- Color coded
- Available in sizes 24 – 14G
- Individually packed

Intraosseous Needle

Product Code: INTRAOS.NDL

- Available in size: 18G, 16G, 15G and 14G
- Suitable for proximal tibial access
- High quality steel stylet which provides excellent bone penetration
- Adjustable depth mechanism
- Buy these needles as singles or as part of the CritiPack® Emergency Intraosseous Pack



[Back to Contents Page](#)